

Colchester Research Facility
Laboratory for Clinical Biochemistry Research (LCBR)
University of Vermont
Department of Pathology and Laboratory Medicine



Lab Director: Russell P. Tracy, Ph.D.


Coordinator: Elaine Cornell

MESA 7 Project Manager: Jessica Rooney, MPH

MESA Repository Manager: Rebekah Boyle, MS

LCBR

- Russell P. Tracy, Ph.D has directed the LCB since its inception in 1986; ABCC Board-certified Clinical Chemist
- Current Personnel includes 8 faculty members, 25 technical and administrative staff, plus students, post-doctoral fellows, and visiting scientists
- Focus on epidemiological & clinical trials research in the areas of coagulation, fibrinolysis, fibrosis, and innate and adaptive immunity
- Frequently serve as a central analysis laboratory and sample repository



Role of the LCBR in MESA 7

We are the MESA 7 Central Blood Analysis Lab (CBAL) and Repository



What we do:

- Assist with planning and development of MOPs for blood collection, handling, storage, sample utilization, etc.
- Receive shipments of participant samples from study sites
- Provide quality assurance feedback
- Perform some biomarker testing along with the University of Minnesota
- Inventory and provide repository storage for study samples

Blood Collection



Set Up – Venipuncture Supplies



Draw Tubes


MESA 7 Draw			
(All Study Sites Except Northwestern)			
Quantity	Type	Volume	Total Volume
4	EDTA	10 mL	40 mL
1	Serum	10 mL	10 mL
1	Serum	5 mL	5 mL
1	EDTA	2 mL	2 mL
1	EDTA	10 mL	10 mL
1	PaxGene RNA	2.5 mL	2.5 mL
	Total Draw		69.5 mL

MESA 7 Draw





Draw Tubes Continued

MESA 7 + Epi Draw (Northwestern)			
Quantity	Type	Volume	Total Volume
4	EDTA	10 mL	40 mL
1	Serum	10 mL	10 mL
1	EDTA	2 mL	2 mL
1	EDTA	6 mL	6 mL
1	CPT	8 mL	8 mL
1	PaxGene RNA	2.5 mL	2.5 mL
	Total Draw		68.5 mL



Aliquot Tube Preparation

Label Orientation



Safety Issues and Precautions for Handling Blood Specimens

In accordance with the OSHA regulations on bloodborne pathogens, the LCBR recommends the following lab safety protocol for the field center laboratories:

- ✓ Use of non-permeable lab coats, latex gloves, and face shields when handling any blood in any situation where splashes, spray, spatter, or droplets of blood may be generated and eye, nose, or mouth contamination can be reasonably anticipated.
- ✓ Use of aerosol containers in all centrifuges.
- ✓ Follow 'Standard Precautions' when handling any blood products.
- ✓ Contaminated needles and sharps shall be immediately placed in a puncture-resistant, leak-proof container. **Never recap or break needles.**
- ✓ Hepatitis B vaccine should be offered to all unvaccinated technicians handling blood, and documentation of vaccination, or technician's declining to be vaccinated, should be kept on file at the Clinical Center.
- ✓ Limit distractions - Noise, radios, phones, etc for your safety and the safety of those around you

Sources of safety information:

<http://www.cdc.gov/niosh/topics/bbp/>

<http://www.osha.gov/SLTC/bloodbornepathogens/index.html>

<http://www.cdc.gov/elcosh/docs/d0300/d000378/d000378.pdf>



Preparation of Venipuncture Site





Venipuncture procedure





Venipuncture procedure continued



Venipuncture continued




← Remove needle and apply pressure
w/ gauze pad

Apply
bandage →




MESA 7 Classic Samples after Draw




- EDTA - After filling, these tubes are mixed for ~30 seconds on the tube rocker then placed on wet ice. EDTA tubes for plasma preparation are centrifuged within 30 minutes of the draw. EDTA for whole blood can remain on ice until it can be stored in -80 C freezer.
- Serum - After filling, gently invert tube 5 times to ensure proper mixing of the activator with the blood. Tubes remain upright at room temperature for a minimum of 60 minutes, but not longer than 90 minutes, to allow the blood to clot.
- Paxgene - After filling, gently invert tube 8-10 times to mix blood with RNA stabilization additive. The tube can remain at room temperature until it is placed in a -80 C freezer for storage.

Blood Processing




MESA Exam 7 Processing Guide

After the visit, you should have:



Prepare for Centrifugation

Balance tubes to be centrifuged



Use sealed cups to minimize aerosol hazard





Run at 4°C for 30,000 G/min




Pool and Gently Mix

- Have labeled 15 mL (serum) 50 mL (plasma) pooling tubes ready
- Draw plasma or serum off using a pipette
- Use care to avoid disturbing the cell layer at the bottom of the tube
- Cap and gently rock tube to mix




Pour Over Packed Red Blood Cells


- Have labeled transfer tubes ready for 5 EDTA draw tubes
- Carefully pour any remaining plasma and all packed cells into transfer tube
- Use one transfer tube per draw tube
- Cap and freeze upright in labeled 5x5x3 inch box



Aliquot and Cap



Example rack from Exam 6




Completed blood and urine set

- Place rack with pre-labeled cryovials on ice or cold pack
- Aliquot 0.5 mL plasma into first 10 cryos ('01-'10)
- Aliquot 0.25 mL plasma into next 12 cryos('11-'22)
- Aliquot 0.5 mL plasma into remaining 32 cryos('23-'54)
- Aliquot 0.5 mL serum into first 14 cryos ('55-'68)
- Cap full cryos with colored cap (EDTA-blue, serum-red, urine-yellow)
- Change pipette tips between sample types and participants

Processing Notes

- Serum sample processing follows the same steps as EDTA plasma
- Serum must sit for a minimum of 60 minutes before centrifugation to allow for blood to clot
- Fill tubes in **numerical** order
- Do not attempt to make all aliquots if only a partial volume was collected, fill cryos in order with specified volume leaving the last with partial volume
- Mark partial volume samples on the label with a “P” and record on the processing form
- Record all aliquots made and any deviations from protocol on processing form





Urine Collection and Processing



'spot' urine collection
~11 mLs required for
processing



1.5 mL x 6 cryos



Labeled urine aliquots

- Cryovials 69-74, 2 mL size vial
- 1.5 mL urine added to each tube
- Place all tubes on ice until the can be frozen @ -80 C
(preferably within 10 mins)

Prepare for Freezer



- Place all cryovials in 5x5x2 inch fiberboard box with 100 cell grid
- Label box with “Freezer Box” label
- Place Paxgene RNA draw tube, EDTA 2 mL whole blood draw tube, and 5 pRBC transfer tubes in 5 x 5 x 3 inch fiberboard box with 49 cell grid
- All samples should be frozen upright in a -80°C freezer until ready to ship

Packing and Shipping



Packing Frozen Blood Samples

The following shipping instructions comply with the International Air Transport Association's Goods Regulation-Packaging Instructions 650 and 954

- Line Styrofoam mailer with cushioning/absorbent material (lab mat)
- Place 5-10 lbs of dry ice in bottom of mailer
- Place another layer of absorbent material on top of dry ice to act as a barrier between the dry ice and freezer boxes
- Ensure cryovials are packed in 5x5x2 inch cardboard freezer boxes with either 81 or 100 cell grids
- Secure cryovial boxes with rubber bands and place in a zip top bag with absorbent material sufficient for the volume of sample being sent
- Ensure tall tubes are packaged according to instructions and placed in zip top bag with absorbent material sufficient for the volume of sample being sent
- Place bagged samples into Styrofoam mailer
- Place another layer of absorbent material over freezer boxes
- Place remaining dry ice (~5-10 lbs) on top of absorbent material and put the top of the Styrofoam in place
- Insert the shipping form and participant processing forms on top of the Styrofoam mailer before sealing the outer cardboard flaps






Shipping Schedule

- Frozen samples should be shipped Monday or Tuesday of every other week
- We would like to develop a shipping schedule so every site ships the same week, but are divided between the two days
- Please let us know if you have any limitation at your site that are related to schedule (Dry Ice availability, etc)



Shipping Frozen Blood Samples



- Use Fed-Ex Priority Overnight Delivery
- Shipment information should be e-mailed to Jessica Rooney, Elaine Cornell, and Rebekah Boyle (Jessica.lanzer@med.uvm.edu, elaine.cornell@med.uvm.edu & Rebekah.boyle@med.uvm.edu)
- All packages must be marked with
 - A UN3373 Biological Substance Category B label
 - A dangerous goods Category 9 UN1845 Dry Ice label (with dry ice weight)
 - A completed Fed-Ex air bill
 - The shipper and recipient's name, address, and phone number
- All packages should be shipped to:
Rebekah Boyle
University of Vermont
Laboratory for Clinical Biochemistry Research
360 South Park Drive
Colchester, VT 05446
p:(802)656-8938

Certification





Certification Process

Obtaining Certification

- Read and understand the appropriate chapters in the MESA 7 Manual of Operations.
- Attend Central Training Zoom session
- Successfully complete the written exam
- Successfully complete blood draw/processing while be observed

Maintaining Certification

- Perform phlebotomy/processing on a minimum of one person every month for the first three month and every other month after that

Certifying Other Technicians

- Ensure new technician has read MESA's Laboratory MOP and section of the Exam 7 MOP.
- Perform phlebotomy/processing with new technician observing
- Supervise as new technician drawing/processing blood samples from one volunteer as described in the protocol and in accordance with the certification checklist
- Administer written examination prepared by LCBR.
- Send exams and checklists to LCBR and CC



Questions

A wide-angle photograph of a sunset over a large body of water, likely a bay or fjord. The sky transitions from a deep blue at the top to a bright orange and yellow near the horizon. The sun is partially obscured by a thin layer of clouds, creating a shimmering reflection on the water's surface. In the foreground, a wooden deck made of horizontal planks leads towards the water. Two dark, rectangular concrete pillars stand on either side of the deck. The overall mood is serene and contemplative.