

Processor ID: _____

Participant ID: _____

Date: ____ / ____ / ____

Processing Start Time: ____:____ EDTA

Processing Start Time: ____:____ Serum

Processing Start Time: ____:____ Urine

(Please check box below if PAXgene RNA draw tube is collected & frozen)

Cryo #	Type	Color	Sample Vol	Check if Done	Comment*	Cryo #	Type	Color	Sample Vol	Check if Done	Comment*
1	EDTA WB	O	0.5 mLs		(HbA1c WB)	22	Urine	Y	1.5 mLs		
2	EDTA	P	0.5 mLs			23	Urine	Y	1.5 mLs		
3	EDTA	P	0.5 mLs			24	Urine	Y	1.5 mLs		
4	EDTA	P	0.5 mLs			25	Urine	Y	1.5 mLs		
5	EDTA	P	0.5 mLs			26	Urine	Y	1.5 mLs		
6	EDTA	P	1.0 mLs			27	Urine	Y	1.5 mLs		
7	EDTA	P	1.0 mLs			28	PAXgene RNA	R	2.5 mLs		(10 mL Draw Tube)
8	EDTA	P	1.0 mLs								
9	EDTA	P	1.0 mLs								
10	EDTA	P	1.0 mLs								
11	EDTA	P	1.0 mLs								
12	Serum	R	0.5 mLs								
13	Serum	R	0.5 mLs								
14	Serum	R	0.5 mLs								
15	Serum	R	0.5 mLs								
16	Serum	R	1.0 mLs								
17	Serum	R	1.0 mLs								
18	Serum	R	1.0 mLs								
19	Serum	R	1.0 mLs								
20	Serum	R	1.0 mLs								
21	Serum	R	1.0 mLs								

* P = partial volume
 O=orange, P=purple, Y=yellow, R=red.

Comments:

LCBR Rec'd Date: _____

Frozen: Y N